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Instructions for use and care Dynair® Balance Master

Congratulations on purchasing your new Dynair® Balance Master.

With the Dynair® Balance Master, you have chosen a TOGU® brand product that is “made in Germany”! It is manufactured in accordance with applicable EU regulations, contains no heavy metals, and is made from 100% recyclable plastic.

Please read the instructions for use carefully.

The Dynair® Balance Master is delivered optimally inflated.


Inflation:

Place the needle on the pump. Hold the material around the valve between your thumb and index finger and carefully insert the slightly moistened needle tip straight into the valve.


Deflating:

Pinch material around the valve with two fingers and moisten the tip of the needle. Insert the needle straight and with caution. Deflate to sired level. (Pump not necessary for this procedure.)


Too much air:

Too little air:

Optimal filling quantity:

Equipment details:

Dimensions: 40 x 40 x 12 cm

Weight: approx. 3840 g

Load capacity: approx. 200 kg

Temperature resistance:

The Dynair® Balance Master is UV-resistant and also heat-resistant to a certain extent. The Dynair® Balance Master can be used without any problems at temperatures up to 45°C. However, the Dynair® Balance Master should not be stored in direct sunlight (e.g., directly on a window pane). The material becomes softer when exposed to heat.

The Dynair® Balance Master is also cold-resistant, but the storage temperature should not fall below -10°C. The material becomes harder when exposed to extreme cold.

Care:

Any mild standard cleaning agents and disinfectants can be used to clean the Dynair® Balance Master. Solvent-based, acidic or caustic cleaners will damage the material.

Avoid contact with sharp or sharp-edged objects that could harm the product.

Use only the pump provided with the product. Other needles might damage the valve.

Dynair® Balance Master adapts to room temperature in only a few minutes.