MD Dynair® Ballkissen® XXL

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TOGU GmbH

Atzinger Str. 1

D-83209 Prien-Bachham

Tel.: +49 (0)8051 9038-0

info@togu.de

Version 2/22

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Medical device class 1

Table of contents

1	General information	3
1.1	Definitions	3
1.2	Abbreviations	5
2	Safety	6
2.1	Safety instructions	6
2.2	Product lifetime	6
2.3	Remaining residual risks	6
2.4	Cautions/Warnings	6
3	Device description	7
3.1	Intended use	7
3.2	Device description	8
3.3	Performance description	8
3.3.1	Product performance	8
3.3.2	Product lifetime	8
4	Application	8
4.1	Commissioning	8
4.2	Correct use	8
4.3	When should the product not be used?.....	9
5	Cleaning and disinfection	9
6	Storage and transport	9
7	Waste disposal	10
8	Legal notice	10

1 General information

Do not use the device until you have read and understood the operating instructions.

All information contained in these instructions for use is the property of **TOGU GmbH** (= manufacturer) and may not be reproduced or duplicated without prior permission. The manufacturer reserves the right to change or completely replace this instruction manual at any time and without prior notice.

You should always ensure that you are working with the latest version of the operating instructions. If there is any doubt about the up-to-dateness, contact the manufacturer.


Keep these instructions for use in a safe place and pass them on to other users of the product before use.

TOGU® and Dynair® Ballkissen® XXL are protected word marks of Gebr. Obermaier oHG.

1.1 Definitions

In these operating instructions, special information is marked as follows:

Symbol	Meaning
	<p>Manufacturer</p> <p>In addition to this symbol, the name and address of the manufacturer, i.e. the person who places the (medical) device on the market, must be indicated</p>
	<p>Date of manufacture</p> <p>Next to this symbol, the date of manufacture of the (medical) product must be indicated in the form YYYY-MM, e.g. 2019-12 for December 2019.</p>
	<p>Usable until</p> <p>For (medical) products that are restricted in their period of use, this symbol shows the date after which the (medical) product may no longer be used in the form YYYY-MM-DD, e.g. 2020-12-01 for 01.12.2020. The date must appear next to the symbol.</p>
	<p>Production lot number/batch</p> <p>Indicates the lot number to which the respective (medical) product belongs. It consists of a 6-digit code and indicates the date of manufacture in the format YYMMDD (year, month, day, each with 2 digits). The lot number is located on the packaging and must be next to the symbol.</p>
	<p>Article number/product identification number</p> <p>Specifies the article or product identification number that must be next to the symbol.</p>

Symbol	Meaning
	Protect from sunlight! If the (medical) product and/or its packaging can be negatively affected in its properties by sunlight, this symbol must be applied to the packaging as a warning.
	Store in a dry place. If the (medical) product and/or its packaging may be negatively affected in its properties by direct moisture, this symbol must be applied to the packaging as a warning.
	Temperature limit. If the (medical) product and/or its packaging can be negatively affected in its properties by certain temperature ranges, this symbol with indication of the lower and upper temperature limits must be affixed to the packaging as a warning. The limit values must be indicated in each case next to the horizontal lower or upper line.
	Air humidity limit. If the (medical) product and/or its packaging can be negatively influenced in its properties by certain air humidity ranges, this symbol with indication of the lower and upper air humidity limits must be affixed to the packaging as a warning. The limit values must be indicated in each case next to the horizontal lower or upper line.
	Air pressure limit. If the (medical) product and/or its packaging can be negatively influenced in its properties by certain air pressure ranges, this symbol with indication of the lower and upper air pressure limits must be affixed to the packaging as a warning. The limit values must be indicated in each case next to the horizontal lower or upper line.
	Observe the instructions for use! If prior reading of the operating instructions is required for correct use of the (medical) product, this symbol must be clearly visible.
	"Reference to eIFU" shall be replaced with the website address where the user can find the electronic user documentation. 
	Attention! If it is absolutely necessary to read the operating instructions before using the (medical) product without risk, this symbol must be clearly visible.
Max. xxx kg	Performance information relevant for the safe use of the product
	MD = Medical Device Marking that the product is a medical device.

Symbol	Meaning
	<p>CE conformity marking (minimum height 5mm). Indicates compliance with the Essential Requirements of the applicable European Directive/Regulation. For Class I medical devices, only the two letters are affixed. For medical devices > Class I, the two letters are followed by the four-digit number of the Notified Body that carried out the conformity assessment procedure.</p>

1.2 Abbreviations

Shortcut	Meaning
JJMMTT	Year (2 digits), Month (2 digits), Day (2 digits)

2 Safety

2.1 Safety instructions

2.2 Product Lifetime

If the service life of 2 years from first use, or at the latest 10 years after manufacture (see lot number) is exceeded, safe use of the product is no longer guaranteed. Replace the product immediately.

2.3 Remaining residual risks

Please keep the exercise field free of sharp objects or edges. Before each use of your product, check that the surface is undamaged. Product safety is no longer guaranteed if the surface is damaged or if the safety limits described in this manual are exceeded.

2.4 Cautions/Warnings

The product must not be misused, modified or changed in any other way. Safe use is no longer guaranteed in this case.

Only use the product on a stable, horizontal, dry and clean surface. In this way, prevent the product from slipping away.

Keep the product away from:

- sharp and pointed objects
- fire or other sources of heat
- Chemicals

Do not wear clothing or shoes with sharp or pointed parts while using the product.

Keep the area around the product free of objects (furniture, fixtures) within a radius of the user's body length. Keep your distance from other persons. This will prevent injuries in the event of a fall.

The product is unstable. You may suddenly lose your balance. Move onto the product carefully and slowly. Use safe, firmly anchored aids to avoid a fall and injuries.

Sweating or moisture from outside (e.g. rain) can make the product slippery. Also check the surface of the product during use. In this way, avoid slipping and injury in the event of a fall.

The product should be used at room temperature. The product can be used in a temperature range from **-10°C to +45°C**. It must not be used in direct sunlight (e.g. by the window) or near any other heat source (e.g. by the heater).

The material is thermoplastic: this means that it becomes softer with increasing heat and firmer with increasing cold.

Check the product for damage, wear, and functionality before use. If there are any deviations, the product must no longer be used. Remove and replace the product if it shows signs of damage. Do not attempt to repair the product.

If the instructions for use are not followed, an MD Dynair® Ballkissen® XXL may burst when struck by a pointed or sharp object.

Seek professional advice (doctor, therapist) before performing exercises with the product. Children, elderly or impaired persons should use the product only under the supervision of an adult, physician or therapist familiar with the product.

In these cases, the product should not be used: see section 4.3.

3 Device description

MD Dynair® Ballkissen® XXL

With its 50 cm diameter, the MD Dynair® Ballkissen® XXL challenges every exerciser to high performance and improves coordination and strength. Thus, you have even more possibilities for a safe sensorimotor and coordinative training. Like all Dynair® Ballkissen® XXL, it is equipped with Senso nubs on one side and a velvety soft surface on the other.

3.1 Intended use

The MD Dynair® Ballkissen® XXL is a medical device and is used as a non-active therapeutic device for:

- Diseases of the musculoskeletal system,
- as prophylaxis in the professional field and
- as prophylaxis in the private/amateur sector.

The product is used by physiotherapists (or similar qualification) in rehabilitation or physiotherapy facilities (or similar), as well as by patients (age from 6 years) in the home environment.

In therapy, the MD Dynair® Ballkissen® XXL is used in orthopedics and neurological rehabilitation, among others, to train mobility, strength and coordination and also serves to promote sensory function, coordination and the musculoskeletal system.

Contraindications :

The MD Dynair® Ballkissen® XXL has been developed for various therapeutic and sporting applications and must not be used for any other purpose than that stated.

The use of the product is prohibited in case of:

- Severe osteoporosis
- Other diseases that greatly reduce the strength of the bony system

Conditional contraindications are defined as follows (only after consultation with physician/therapist):

- Balance disorders
- Dizziness
- Epilepsy

Do not train with weights or dumbbells of any kind on the MD Dynair® Ballkissen® XXL

3.2 Device description

Description	Size	Weight net
MD Dynair® Ballkissen® XXL	50 cm	2500 g

The package includes packaging, quick start guide, instructions with examples of use.

3.3 Performance description

In the table below, the load capacity is specified in relation to the maximum load capacity, i.e. how much weight an MD Dynair® Ballkissen® XXL can safely support in undamaged condition (max. load).

3.3.1 Product performance

Description	max.load
MD Dynair® Ballkissen® XXL	ca. 200 kg

3.3.2 Product lifetime

If the service life of the product exceeds 2 years from the date of first use, or a maximum of 10 years from the date of manufacture (see lot number), safe use of the product is no longer guaranteed. Replace the product immediately.

4 Application

4.1 Commissioning

When opening the packaging, take care not to damage the product.

Check the product for damage, wear and functionality before use. If there are any deviations, the product must no longer be used. Remove and replace the product if it shows any damage. Do not attempt to repair the product.

4.2 Correct use

Observe the detailed safety instructions under point 2.

The product should be used at room temperature. The product can be used in a temperature range from **-10°C to +45°C**. Keep the area around the product free of objects (furniture, fixtures, etc.) within a radius of the user's body length. Keep a distance from other persons. This will prevent injury in the event of a fall.

Before each use, check the product for damage, wear and tear and proper functioning. If there are any deviations, the product must no longer be used. Remove and replace the product if it shows signs of damage.

If the service life of the product exceeds 2 years from the date of first use, or a maximum of 10 years from the date of manufacture (see lot number), safe use of the product is no longer guaranteed. Replace the product immediately.

For the application, please note:

- Adopt a correct posture
- Perform slow and controlled movements
- If pain occurs, seek advice from a doctor or therapist
- Children should only exercise under adult supervision

Inflate:

The MD Dynair® Ballkissen® XXL is delivered optimally filled with air. You should always make sure that the surface is "flat". If you do want to add some air, please use the enclosed hand pump with the black plug valve attachment.

4.3 When must the product not be used?

Contraindications:

The use of the product is prohibited in case of:

- Severe osteoporosis
- Other diseases that greatly reduce the strength of the bony system

Conditional contraindications are defined as follows (only after consultation with physician/therapist):

- Balance disorders
- Dizziness
- Epilepsy.

Other restrictions (expiry of service life see point 3.3.2)

5 Cleaning and disinfection

Clean the product after each use exclusively with solvent-free, non-acidic and non-caustic cleaning agents and disinfectants. The device is waterproof.

6 Storage and transport

Store the product in a place free from sharp and pointed objects to avoid damaging the product.

Da material is thermoplastic, which means it softens with increasing heat and solidifies with increasing cold. The product should be stored at room temperature. The product can be stored in a temperature range from -10°C to +45°C.

It must not be stored in direct sunlight (e.g. at the window) or at any other heat source (e.g. at the heater).

Use the supplied packaging when transporting the product.

7 Waste disposal

Products that are no longer used can be returned freight prepaid to

TOGU GmbH
Atzinger Str. 1
D-83209 Prien-Bachham

be returned and are fed into the recycling process.

8 Legal notice

The manufacturer accepts no responsibility for personal injury, damage to animals, property or the environment resulting from improper use. Please always follow the instructions and warnings accompanying the product.

Report serious incidents involving the use of the product to the manufacturer and the competent authorities of the Member State.

On the product packaging there is a 6-digit code preceded by the abbreviation LOT. This code indicates the date of manufacture in the format YYMMDD.

In the event of a complaint, this code must be used to identify the batch and the materials used.

Classification: Medical device class 1