Übungsposter Senso® Balance Igel

Exercise Poster Senso® Balance Hedgehog

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GB Lying on one side, lower leg flexed at right angle, upper leg extended. A Senso® Balance Hedgehog is placed under the supporting arm, i. e. under the elbow. Support body with arm and raise the pelvis, then lower pelvis again until it contributes to support.

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Lying on one side, placed a Senso® Balance Hedgehog under the extended lower arm and another Senso® Balance Hedgehog under the flexed knee. The foot touches the floor with the upper leg extended. Hold upper arm extended above head. Pull body together laterally: raise the upper leg, at the same time raising the trunk, whereby the upper arm with the hand extended upright moves in the direction of the heel. Then return to first position.

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On all fours on the Senso® Balance Hedgehogs (Senso® Balance Hedgehogs either under hands and/or knees or diagonally):
- extend right or left arm
- extend right or left leg
- extend right arm and left leg and vice-versa

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On all fours on the Senso® Balance Hedgehogs, whereby the Senso® Balance Hedgehogs under the hands are a bit farther apart and the feet are folded over. Lower the trunk and at the same time flex the elbows, remaining within the arm axis, then return to first position by re-extending the arms.

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Set up the Senso® Balance Hedgehogs in a row, a circuit or randomly. Walk on the Senso® Balance Hedgehogs forwards, backwards and sideways.

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Stand on a Senso® Balance Hedgehog, lift one leg and write with it numbers, letters or words in the air.

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Sit on a Senso® Balance Hedgehog, lift both legs and balance, at the same time
- passing a ball over head
- passing a ball under one leg or both legs
- throwing up and catching a ball

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
Sit on Senso® Balance Hedgehog with the hands supporting you on the floor behind your back. Lift alternately the right leg, the left leg, the right hand, then the left hand. Also try this without the supporting hands.

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On all fours, place a Senso® Balance Hedgehog on the back, then
- crawl across the room
- extend one arm and/or one leg
- crawl towards one another, shake hands, all without losing the Senso® Balance Hedgehog.

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Line up Senso® Balance hedgehogs in two rows. Stand behind one another on the Senso® Balance Hedgehogs and pass a ball or Senso® Balance Hedgehog from one to the other:
- above the head;
- through your legs;
- around to the side.

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Place Senso® Balance Hedgehogs in a circle. Stand within the circle and:
- throw; - bounce; - or pass balls around to the right and left around the circle, whereby you occasionally run around to a new position on the Senso® Balance Hedgehogs; - dribble the ball while walking around on the circle of Senso® Balance Hedgehogs.

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Sit on the ground with Senso® Balance Hedgehogs under the feet, under the hands and possibly also under the buttocks. Raise and hold the body on feet and hands, then lower again.

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Exercise Poster Erizo Senso® Balance

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© TOGU® GmbH, Atzingerstr. 1, 83209 Prien - Germany