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
Manual and Maintenance Instructions Aero-Step® Pro actisan

Congratulations on your purchase of the Aero-Step® Pro actisan.

With the Aero-Step® Pro actisan you have chosen a TOGU brand product made in Germany.

It was produced as a sports device in compliance with valid EU norms, it contains no heavy metals and is made of 100% recyclable plastic.

Please read carefully through the following tips and all aspects of the Aero-Step Pro actisan so that you will be able to enjoy using the product for as long as possible.

The new Aero-Step is more dynamic, bigger and versatile. You can now use both sides – 2 in 1, ideal to clean and store through easy store and easy clean.

The approved tool, with the twin air chamber system provides training with maximal effect. More intensive training for the deep lying musculature apperception, processing of stimuli and stability will be trained and enhanced. Applicable. Very variable, applicable using on standing, kneeling, or lying down very robust. For using in rehabilitation, fitness or at home. Over 1000 soft knobs massage to enhance blood circulation.

How actisan works?

The secret: Actisan is an ingredient for the synthetic material used for the tools and highly effective against bacteria and other microorganisms on an ionic basis. It is constantly building up ions which are harmless for humans. These ions are attacking only the metabolism of the harmful organisms. What is remarkable about Actisan: It is antibacterial as well as antimicrobial and destroys bacteria and other microorganisms at the same time. Therefore, it is much more efficient than the usual wipe off disinfection and this effect lasts for at least five years.

Equipment details:

Your Aero-Step Pro actisan is manufactured from high quality RUTON. This guarantees that the product will be long lasting. Please do not try to fold the Aero-Step Pro actisan.

Maintenance of Aero-Step Pro actisan:

For cleaning you may use any of the mild cleaning and disinfecting substances to be found in most shops. Cleaning substances containing solvents are damaging to the material of which the Aero-Step Pro actisan is made. The Aero-Step Pro actisan is waterproof. You should, nevertheless, ensure that the Aero-Step Pro actisan is dry before you use it, since it may slip if the surface on which it rests is wet. It is not necessary to disinfect the Aero Step Pro from the actisan Studio Line.

Temperature-Resistance:

Aero-Step Pro actisan is resistant to ultra-violet light and resistant to heat, within certain limits. The Aero-Step Pro actisan can be used problem free in temperatures up to 45°C. It should be noted, however, that the higher the temperature becomes the softer the Aero-Step Pro actisan material becomes. In summer temperatures you should not place the Aero-Step Pro actisan directly on an asphalt surface because asphalt stores the heat so much (this is particularly important in Mediterranean countries, in places such as leisure facilities). The Aero-Step Pro actisan is also resistant to cold. It is quiet possible to store away the Aero-Step Pro actisan at temperatures of down to -10°C. It should be noted, however, that the colder the temperature becomes the harder the material of the Aero-Step Pro actisan becomes.

Using the Aero-Step Pro actisan:


The Aero-Step Pro actisan can be used on almost any kind of surface. The patented design of the Aero-Step Pro actisan ensures a sure and safe surface contact. You should, however avoid using the Aero-Step Pro actisan on stony ground (gravelled, pebbled) or on asphalt. Although it is waterproof the Aero-Step Pro actisan is, nevertheless, unsuitable for use in water because it floats. When using the product ensure there are no sharp objects or sharp edges in the immediate vicinity because these can damage the Aero-Step Pro actisan.

Weight limit for the Aero-Step Pro actisan:

The Aero-Step Pro actisan can carry weights of up to 200 kg. It should be noted, however, that if the Aero-Step Pro actisan is being subjected to such upper weight-levels then the air pressure should be adjusted more often because otherwise the material can become fatigued more quickly.

Inflation tips for the Aero-Step Pro actisan:

The Aero-Step Pro actisan has a 2-chamber system. Each of the two chambers is equipped with a needle valve. After longer periods of use it is possible to re-adjust the air pressure yourself. To do this, please use the inflating needle of the TOGU® ball pump. This needle can be used to inflate and to deflate. You can also use any basketball or volleyball needle without back-pressure valve. Moisten the needle before use and make sure that you introduce the needle straight into the valve. To ensure the continued flexibility of the Aero-Step Pro actisan you must take extra care that you do not pump in too much air. The stepping-surface of the Aero-Step Pro actisan should not bulge outwards because if it does so the effectiveness of the training will be greatly reduced (see sketch). Furthermore, if the surface bulges then this means that the Aero-Step Pro actisan is no longer lying flat on the ground or floor, and the stabilising plate may work free of the Aero-Step Pro actisan. Please make sure that both chambers are filled equally, so that an ideal training result can be reached. When used a lot, the Aero-Step Pro actisan may loose some air pressure. In that case just inflate it again (like described above).

Exercises available at: www.togu.de