**MFT Bodyteamwork App oder Software zum Download for free unter: [www.mft-bodyteamwork.com/support](http://www.mft-bodyteamwork.com/support) und AppStore, Google Play, Amazon Appstore**

**MFT Bodyteamwork GmbH,**  
Binderfeld 6, A-6365 Kirchberg/Tirol  
office@mft-bodyteamwork.com  
[www.mft-bodyteamwork.com](http://www.mft-bodyteamwork.com)

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## MFT Trainings Devices Instruction Manual

**Intended Usage:** MFT training devices can be used in athletic and fitness training to improve ones coordination and proprioception and in Physical Therapy for injury prevention and rehabilitation of ones ankles, knees, hips and spine.

**Requirements:** Constant supervision and instruction from a doctor, certified trainer or therapist is required when using MFT training devices for Physical Therapy.  
**Age:** from 4 to 100 years, according to the level of fitness.  
**Weight:** max. 120 kg body weight (approx. 270 lbs)  
**Skills:** no special skills necessary

**Safety instructions:** One must execute all workouts on a compact, non slippery, dry and level surface must be clean and dry during all training. The product cannot be combined with other products. This product must not be used as a toy. Changes to the device are not permitted. Ensure that there is no free fall space and, if necessary, place fall-resistant underlays on top. Children may only train under the supervision of an adult person.

**Additional information:** The device is not suitable for jumping or standing, but only for slow, sensitive movements while sitting.

**Regulations/Maintenance:** Before usage check the device for damages.

DIN EN 50419  
According to 2001/95/EG  
WEEE-Reg.Nr. DE65255935 (ElektroG)  
Reg.Nr.21000465 (BattG)  
UN-Nr. 3091  
<https://www.mft-bodyteamwork.com/CE/>

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Google Chrome, Amazon Fire TV, etc.

## Important information:

**Software-download / more information:**

**[www.mft-bodyteamwork.com/en/support](http://www.mft-bodyteamwork.com/en/support)**

**App download:** AppStore, GooglePlay, Amazon Appstore

### Commissioning

• Remove the battery fuse strip before first use (photo 1)

### Activation of MFT Balance Sensor with Bluetooth

• Activate Bluetooth on the device (tablet, smartphone, computer)

• Move the cushion slightly – tap – start app/software.

• If multiple devices are detected, select the appropriate device.

### Calibration & seating position

• The cushion must be on the chair (bench etc.) when starting the app

• Sit on the middle of the cushion with both feet on the ground

• The balance sensor is correctly positioned when the MFT logo is

pointing forward – between the legs and towards the terminal (see photo 2)

• During training, the location of the cushion should not be changed,

otherwise it must be recalibrated.

### System requirements:

Bluetooth from 4.0 LE, Apple from iOS 6, Mac OSX from 10.7, Android from 4.3,

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