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Usage and care instructions for the Backswing® Trainer

Congratulations on purchasing your new Backswing® Trainer.

With the Backswing® Trainer, you have chosen a TOGU® branded product that is 'made in Germany'!

As a piece of sports equipment, it is manufactured in line with the applicable EU regulations, it contains no heavy metals and is made of 100% recyclable plastic.

The Backswing Trainer is a versatile training device and represents a further development of the proven Airgo lordosis back cushions. An even and secure connection to the base is achieved through the additional stable plate, to which the air cushion is attached. The device can be used on the floor while standing, kneeling, and in further positions to train balance, coordination and stability. Furthermore, the device can be used as a base for patients on their stomach, side and back for mobilisation. Also ideal as an unstable base for stomach and back training in all target groups. The integrated needle valve makes it possible to adapt the inflation to suit the patient or person working out. Fields of application: fitness, functional training, competitive sports

Warnings:

- Over inflation causes the material to over-expand. In this case, functionality and safety can no longer be guaranteed.
- Avoid exposure to direct sunlight and storage close to heat sources, as the material of the Backswing Trainer can otherwise become deformed.
- Avoid contact with sharp and pointed objects (e.g. glass shards/sharp stones in the soles of shoes), as these can damage the device. Do not throw device!

Brief device information:

The Backswing Trainer is made from high-quality Ruton with a needle valve.

Dimensions: approx. 46 x 29 cm

Load-bearing capacity: approx. 150 kg

Maintenance:

To clean the Backswing Trainer, you can use all mild, commercially-available cleaning agents and disinfectants. Cleaners containing solvents, acids or corrosive substances attack the material. The Backswing Trainer is waterproof. Pay particular attention that the underside, which gives the Backswing Trainer its special stability, is always free from dirt and impurities, and is dry.

Inflation:

If desired, the inflation of the Backswing Trainer can be regulated with a commercially-available needle pump (e.g. football pump). The Backswing Trainer is fitted with a needle valve for this purpose. Only make changes to the internal pressure of the Backswing Trainer (inflation or deflation) when the Backswing Trainer is warmed to room temperature.

Over inflation causes the material to over-expand. In this case, functionality and safety can no longer be guaranteed. Always moisten the needle of the pump before carefully inserting it into the needle valve.

Storage:

Please store the Backswing Trainer in a flat, horizontal position.