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User Manual and Care Instructions TOGU Rückenfreund

Congratulations on purchasing your TOGU® Rückenfreund – Made in Germany!

With the Rückenfreund, you have chosen a high-quality and durable training device designed to support **spinal health**, promote **gentle movement**, and help relieve **tension** – sustainably developed and produced in Prien am Chiemsee, Germany.

What is the TOGU Rückenfreund?

A dynamically air-filled tool designed to gently stretch and mobilize the back, torso, and hip area. It specifically targets shortened muscles and fascia caused by prolonged sitting or everyday posture habits – ideal for anyone looking for an easy and time-efficient way to care for their back.

What are the benefits?

- Promotes spinal health and posture
Helps to release fascial adhesions and muscle tension
Improves spinal mobility
Activates deep back muscles through dynamic air movement
Trains the sensorimotor system
Can be used on the upper or lower back
Suitable for use at home, in the office, or while traveling

How does it work?

The Rückenfreund comes pre-inflated and is ready to use.

Simply lie in a comfortable position (refer to the exercise guide) for at least 30 to 60 seconds per exercise. The ergonomic shape, air filling, and high-quality material create gentle, dynamic micro-movements, which are transferred to the body.

These movements stimulate and strengthen deep muscles – a benefit that cannot be achieved with rigid surfaces.

You can adjust the intensity of the stretch using the included pump continuously and individually.

Technical Data

- Dimensions: approx. 30 cm diameter, approx. 10 cm height
Weight: approx. 600 g
Load capacity: up to 300 kg
Material: long-lasting, 100% recyclable plastic
Production: Made in Germany

Care Instructions

- Clean the Rückenfreund with a soft, damp cloth.
Mild, commercially available cleaning agents can be used if needed.
Do not use acidic, solvent-based, or abrasive cleaners.
Store in a clean, dry place at room temperature (storage range: -10 °C to +45 °C).
Avoid direct sunlight and heat exposure.

Important Note

If you have any health issues or experience pain, please consult a doctor or physiotherapist before use.